



GPCR signaling




G protein-mediated signaling by GPCRs





Arrestin-mediated signaling by GPCRs




GPCR biased signaling




GPCR classification



THP-1 cells
CXCR2 antibody [HL2604]
(GTX639056)



SW480 cells
WT CCR4 KO
CCR4 antibody [HL2481]
(GTX638830)




Mouse striatum
mGluR5 antibody [HL3509]
(GTX641397)




Smoothed homolog antibody
[HL2975] (GTX640362)



Dopamine Receptor D2 antibody
[HL1478] (GTX636952)




Mouse testis
FSH Receptor antibody [HL2926]
(GTX640269)




Raji xenograft
[HL2189]
CXCR7 antibody [HL2189]
(GTX638193)


- References:
1. Bone Res. 2025 May 12;13(1):52.
 2. Mechanobiol Med. 2024 May 14;23(1):00078.
 3. Nat Rev Nephrol. 2024 Nov;20(11):722-741.
 4. Br J Pharmacol. 2024 Jun;191(14):2159-2181.
 5. Mol Biomed. 2023 Dec 4;4(1):48.
 6. Pharmacol Rev. 2021 Jan;73(1):89-119.
 7. Nat Rev Mol Cell Biol. 2018 Oct;19(10):638-653.
 8. Nat Rev Mol Cell Biol. 2009 Dec;10(12):919-930.




Human kidney
[HL1864]
GPCR5A antibody [HL1864]
(GTX637589)




CD97 antibody [HL1925]
(GTX637674)




Fritzzled 4 antibody [HL2824]
(GTX640095)




Rat tissues
GABA B Receptor 1 antibody
[GTR3] (GTX641938)



[HL1472]
MC1 Receptor antibody [HL1471]
(GTX636945)



CCR1 antibody [HL3023]
(GTX640436)



HL-60 cells
Gαs Receptor antibody [HL2831]
(GTX640102)

